## **ASOPS**


## Asynchronous Optical Sampling System



In time-resolved measurements an ultrafast pulse triggers a reaction and a second pulse takes a snapshot of the induced change. By shifting the arrival time of the probe pulse with respect to the pump pulse the stimulated process can be followed in time. The asynchronous optical sampling technique allows high-speed scanning over some nanoseconds of time delay without a mechanical delay line. The ultrafast lasers delivering the pump and probe pulses are locked together at a tunable repetition rate difference.

The lasers can be also locked to the same repetition rate value and by shifting the relative phase between the laser pulses, the system will allow measurements in a reduced time window of some 100 ps. Switching between the two modes of operation can be done at the touch of a button.

#### **ASOPS LASER HEAD**



## **MenloSystems**

#### **KEY SPECIFICATIONS**

- Wavelength TWIN: 1560 nm

  DUAL COLOR 780 nm or 1560 nm
- Time Measurement Window 4 ns or 10 ns
- Repetition Rate 250 MHz or 100 MHz

#### **APPLICATIONS**

- Two-Color Pump-Probe Spectroscopy
- Time-Domain THz Spectroscopy
- Material Characterization

#### **FEATURES**

- Faster Data Aquisition Times
- Increased Temporal Measurement Windows
- No Moving Mechanical Components (Improved Beam Pointing Stability, Increased Scanning Speed)
- ASOPS Control Software
   Full Control of the ASOPS Electronics
   via Graphical User Interface, XML-RPC
   Interface Enabling Remote Control,
   Includes PC and Counter

#### OPTIONS

■ VARIO
User-Defined Repetition Rate
Factory-set value selectable in the
50-250 MHz range

#### MULTIBRANCH Additional Seed Ports

Seeding of multiple amplifiers with optional subsequent frequency conversion to cover multiple wavelengths

#### CUSTOM WAVELENGTHS

- Yb Systems (1040nm, 520nm)
- Er Systems (1560nm, 780nm)
- Hybrid Systems

# **ASOPS**



## Asynchronous Optical Sampling System

<150 fs

SPECIFICATIONS	ASOPS TWIN	ASOPS DUAL COLOR
Repetition Rate	250 MHz	100 MHz
Repetition Rate Offset Tuning Range	$\Delta$ f = -10 kHz +10 kHz, in steps of 10 <sup>-5</sup> Hz	$\Delta\text{f} = \text{-10 kHz}$ +10 kHz, in steps of 10-5 Hz
Time Measurement Window	4 ns	10 ns
Scan Duration 1 / $\Delta$ f *	0.1 ms @ 10 kHz offset, 1 s @ 1 Hz offset	0.1 ms @ 10 kHz offset, 1 s @ 1 Hz offset
Data Point Increment **	160 fs @ 10 kHz, 0.016 fs @ 1 Hz	1 ps @ 10 kHz, 0.1 fs @ 1 Hz

#### **SPECIFICATIONS LASERHEADS**

RMS Timing Jitter [0.1 Hz - 500 kHz]

Wavelength	1560 nm	1560 nm	780 nm
Average Output Power	>75 mW (from each laser)	>100 mW	>100 mW
Output Port	fiber-coupled FC/APC	free space	free space
Pulse Length	<150 fs after 6 m PM fiber	<90 fs	<120 fs
Piezo Tuning Range	>625 Hz	>100 Hz	>100 Hz
Piezo Bandwidth	>30 kHz	>30 kHz	>30 kHz
Stepper Motor Tuning Range	>2 MHz	>330 kHz	>330 kHz
Trigger Signal	TTL level at offset frequency, <10 ns rise time	TTL level at offset frequency	uency, <10 ns rise time

<sup>\*</sup>Scales inversely with the repetition rate offset. \*\*Scales with the ratio of the repetition rate offset and the repetition rate squared (\Delta f/f<sub>r</sub><sup>2</sup>).

#### **REQUIREMENTS**

Operating Voltage	110/115/230 VAC	110/115/230 VAC
Frequency	50 to 60 Hz	50 to 60 Hz
Cooling Requirements	no water cooling required	no water cooling required
Operating Temperature	22 ± 5 °C	22 ± 5 °C
Optical Unit Dimensions/Weight	415 x 400 x 110 mm <sup>3</sup> , 35 kg	500 x 535 x 110 mm <sup>3</sup> , 35 kg
Control Electronics Dimensions/Weight	mounted in a 19" rack cabinet, 800 x 600 x 1800 mm <sup>3</sup> , 75 kg	mounted in a 19" rack cabinet, 800 x 600 x 1800 mm <sup>3</sup> , 75 kg

ORDERING INFORMATION		
Product Code	ASOPS TWIN	ASOPS Dual Color

Please call for pricing. Specifications are subject to change without notice. Custom modifications are available, please inquire.

<150 fs





### **MenioSystems**

Menlo Systems GmbH T+49 89 189 166 0 sales@menlosystems.com Menlo Systems, Inc. T+1 973 300 4490 ussales@menlosystems.com **Thorlabs, Inc.** T+1 973 579 7227 sales@thorlabs.com

